



MODEL 190

Cryo-Can

Eliminates sample contamination during scanning electron microscope (SEM) operation.

Model 190 Cryo-Can Specifications

Capacity	270 ml liquid nitrogen
Hold time	Up to 6 hours
SEM operating conditions	Up to 30 keV
Flange interface	Fits all SEMs; subject to the availability of ports
Liquid Nitrogen vessel	Can be baked to 150 °C
SEM compatibility	Contact a Fischione Instruments representative to discuss SEM compatibility
Warranty	One year

Distribution in the UK & Ireland



Lambda Photometrics Limited

Lambda House Batford Mill
Harpenden Herts AL5 5BZ
United Kingdom

E: info@lambdaphoto.co.uk

W: www.lambdaphoto.co.uk

T: +44 (0)1582 764334

F: +44 (0)1582 712084

